



Test Menu

Genetic Assays, Inc.
 4711 Trousdale Drive, Suite 209
 Nashville, TN 37220
WWW.GENETICASSAYS.COM
 PH: (615) 781-0709
 TF: (800) 390-5280
 FX: (615) 781-0766

Code	Assay*	Targets	Methodology	Results	Reference Range	TAT	Specimens*
3701	Adenovirus	DNA	qPCR	Qual/Quant	100 to 1.0 x 10 ¹⁰ copies/mL	same day	Blood, Respiratory, Plasma, Marrow, CSF, Serum, Swab, Stool, Urine
824	Bacterial Vaginosis Panel	7 causative pathogens	rPCR	Qual/ Semi-quant	Not Detected <i>Gardnerella Vaginalis</i> - Normal	24-48 hrs	G Swab®, ThinPrep
3700	BK Virus	DNA	qPCR	Qual/Quant	500 to 1.0 x 10 ¹⁰ copies/mL	same day	Blood & Urine (serial testing) Plasma, Serum, CSF, Marrow
3333	CT (<i>Chlamydia trachomatis</i>)	Ribosomal RNA	Aptima Combo 2® TMA	Qual	Not Detected	24-48 hrs	G Swab®/Urine, ThinPrep, SurePath, Swab, Urine
0180	NG (<i>Neisseria gonorrhoeae</i>)						
6262	Cystic Fibrosis Mutation Detection	39 most common mutations & 4 associated polymorphisms	xTAG™ CF 39 kit v2	Qual/Geno	Interpretive Report	1x/week	Blood, ThinPrep, SurePath, Buccal Swabs
3702	Cytomegalovirus (CMV)	DNA	qPCR	Qual/Quant	200 to 1.0 x 10 ¹⁰ copies/mL	same day	Urine, CSF, Swab (saliva), Blood, Plasma, Fluid, Respiratory, Stool
4025	Enterovirus	RNA	RT-PCR	Qual	Not Detected	same day	CSF, Plasma, Blood, Urine, Serum, Amniotic Fluid, Swab, Stool
6111	Epstein-Barr Virus (EBV)	DNA	qPCR	Qual/Quant	200 to 5,000,000 copies/mL	same day	Blood, CSF, Plasma, Serum, Marrow
7667	Gastrointestinal Pathogen Panel (GPP) (22 Targets)	Bacteria, Diarrheagenic <i>E. coli/Shigella</i> , Parasites, Viruses	FilmArray® multiplex RT-PCR	Qual	Not Detected	1-6 hrs	G Swab®, Raw Stool
219TQ	HCV Viral Load Monitoring	RNA	RT-qPCR	Qual/Quant	12 to 1.0 x 10 ⁷ RNA IU/mL	24-72 hrs	Serum, Plasma
8698	HCV Genotyping	HCV types 1a, 1b, 1, 2, 3, 4, & 5	RT-rPCR	Genotype	Genotype classifications will be reported	1x/week	
875TQ	HIV-1 Viral Load Monitoring	RNA	RT-qPCR	Qual/Quant	40 to 1.0 x 10 ⁷ RNA IU/mL	1-4 days	Plasma, CSF
395H	HPV High Risk w/ ID of 16 & 18/45	DNA of 14 HR types + 16 & 18/45	Hologic® Aptima®	Qual	Not Detected	MWF	ThinPrep, SurePath, Biopsy
478	HPV w/ ID of 16, 18, & 45	HPV DNA (all types) + 16, 18, & 45	PCR + rPCR			1x/week	
7575	HPV Genotyping	HPV DNA (all types)	PCR & DNA Sequencing			1x/week	
900	HSV-1&2 (2 Targets)	HSV type 1 DNA HSV type 2 DNA	rPCR	Qual	Not Detected	same day	CSF, G Swab®, ThinPrep, SurePath, Swab, Blood, Tissue
273	Meningitis Encephalitis Panel (MEP) (14 Targets)	6 Bacteria, 7 Viruses, 1 Yeast	FilmArray® multiplex RT-PCR	Qual	Not Detected	1-6 hrs	CSF
250	Mycobacteria	Mycobacteria DNA (all species) <i>M. tuberculosis</i> complex	rPCR	Qual	Not Detected	MWF	Bronchial Washings, Sputum, Positive AFB Media (Concentrate), CSF, Bodily Fluid, Tissue (Fresh or Paraffin Embedded)
275	w/ AFB Smear & Culture					same day	
250FFPE	Formalin-Fixed Paraffin-Embedded Tissue	All other atypical species	DNA Sequencing	Qual	Species identification	1x/week	
1000	Mycobacteria DNA Sequencing						
2001	Respiratory Virus & Bacteria Panel (RVBP) 20 Targets	17 Viruses 3 Bacteria	FilmArray® multiplex RT-PCR	Qual	Not Detected	1-6 hrs	Upper Respiratory Swab, Bronchial Wash or Lavage
301	STD3 Panel (3 Targets)	CT/NG and <i>Trichomonas vaginalis</i>	rPCR	Qual	Not Detected	24-48 hrs	ThinPrep, SurePath, G Swab®/Urine
501	STD5 Panel (5 Targets)	CT/NG, HSV-1&2, <i>Trichomonas vaginalis</i>	rPCR	Qual	Not Detected	24-48 hrs	ThinPrep, SurePath, G Swab®
8425	Tick-Borne Ehrlichiosis Panel (4 Targets)	1 <i>Anaplasma</i> and 3 <i>Ehrlichia</i> species	rPCR	Qual	Not Detected	same day	Whole Blood, CSF
8667	Tonsillitis Panel (5 Targets)	Adeno, CMV, EBV, HSV-1&2	rPCR	Qual	Not Detected	48 hrs	Throat Swab (G Swab®)
101	<i>Trichomonas vaginalis</i>	<i>Trichomonas vaginalis</i> DNA	rPCR	Qual	Not Detected	24-48 hrs	ThinPrep, SurePath, G Swab®/Urine

* Please contact GA or visit website for detailed test information sheets.

PCR = Polymerase Chain Reaction
 q = quantitative real-time
 r = real-time
 RT = Reverse Transcriptase

* Please contact GA for questions about other acceptable specimens.